







AcuClot Biopolymer Based Haemostatic Sponge

AcuClot is an advanced hemostatic sponge, designed to stop Surgical/ Superficial bleeding. Based on ActClot Technology (Activated Clotting Mechanism), the composition and structure accelerates blood coagulation cascade, resulting in faster bleeding control.



Charge-based platelet aggregation, enhances blood coagulation (Ref. MOA).

- Due to its cationic nature, the sponge is bioadhesive and creates a mechanical seal at the wound and accelerates the natural process of blood clotting (Ref CH T: 1.2, 2.1).
- Microporous honeycomb structure, supports in higher (up to 40 times of sponge weight) and rapid blood absorption (Ref CH T: 1.1, 1.2).
- Sponge is based on Biopolymer*, hence completely biodegradable



Acuclot hemostatic sponge is recommended for being used on:

- Heparinized Patients
- Friable Tissues
- Haemophilia Patients
- Trauma Injuries
- Operative Procedures
- Dental Extractions
- Dialysis Puncture Sites
- Catheter Sites







EASY Every Second Counts



Mechanism of Action:



ActClot Technology: AcuClot plays a Haemostatic role by causing **Platelet Aggregation and Inhibiting Fibrinolysis**

References-

- 1. Bano, Ijaz; Arshad, Muhammad; Yasin, Tariq; Ghauri, Muhammad Afzal; Younus, Muhammad (2017). Chitosan: A potential biopolymer for wound management. International Journal of Biological Macromolecules, 102, 380–383.
 Hu, Zhang; Zhang, Dong-Ying; Lu, Si-Tong; Li, Pu-Wang; Li, Si-Dong (2018). Chitosan-Based Composite Materials for Prospective Hemostatic
- Applications. Marine Drugs, 16(8), 273. 3. Chih-Hsin Wang;Juin-Hong Cherng;Chuan-Chieh Liu;Tong-Jing Fang;Zhi-Jie Hong;Shu-Jen Chang;Gang-Yi Fan;Sheng-Der Hsu; (2021).
- Procoagulant and Antimicrobial Effects of Chitosan in Wound Healing . International Journal of Molecular Sciences.
- 4. Smith, A. H.; Laird, C.; Porter, K.; Bloch, M. (2013). Haemostatic dressings in prehospital care. Emergency Medicine Journal, 30(10), 784-789.
- 5. Khoshmohabat, Hadi; Paydar, Shahram; Kazemi, Hossein Mohammad; Dalfardi, Behnam (2016). Overview of Agents Used for Emergency Hemostasis. Trauma Monthly, 21(1).



Characterization Tests:

AcuClot Haemostatic Sponge

1. Microscopic Evaluation (SEM Analysis):

- Scanning Electron Microscopy image showing highly porous matrix with a typical honeycomb microstructure (100-300 μm).
- 1.2 Uniform porosity allows it to **rapidly absorb** and block plasma during haemostasis.



2. Bioadhesiveness:

- 2.1 AcuClot's cationic nature enables it to adhere strongly to the surface.
- 2.2 Strong bioadhesive forces create a **mechanical seal** accelerating the natural mechanism of blood clotting.



3. Absorbency:

3.1 AcuClot's ability to absorb liquid is up to 40 times its weight and its ability to absorb blood is about 15-20 times its weight.





Comparative Tests:

AcuClot Vs Cotton Gauze

1. Whole Blood Clotting Time Study:

- 1.1 AcuClot quickly forms the blood clot and remains **adhered** to the test tube, unlike cotton gauze.
- **1.2** The formation of a blood clot is confirmed by **inverting** the tubes.



2. Haemolysis:

- 2.1 Blood cells that have not clotted are haemolyzed in water, giving them a characteristic **red colour**.
- 2.2 AcuClot does not show haemolysis as it forms a very strong clot showing it's strength.



How to Apply?



Open the pack and take out **AcuClot.**



Apply AcuClot to the bleeding site so that it completely covers it. Apply consistent pressure until AcuClot adheres to it.



If necessary, use a secondary dressing to keep **AcuClot** in place.



Removal: Irrigate AcuClot thoroughly with water or saline & gently pull it away.

Post Catheterization

Catheter Sites, Dialysis Sites and Any Punctures



Variants: ACT0303 , ACT0505

Dental Procedures

Oral Surgeries, Extraction Sockets and Multiple Tooth Extractions





Variants: ACT0101

Operative Procedures

Surgical Incisions, Resections and Arterial Bleeds



Variants: ACT0104, ACT0106, ACT0108, ACT0505, ACT0808

Trauma Injuries

Lacerations and Abrasions



Variants: ACT1010

	APPLICATION	PRODUCT CODE	SIZE(cm)	UNITS PER BOX
Ordering Information	Dental	ACT0101	1x1	20
	Haemorrhoidectomy	ACT0104	1x4	20
	Haemorrhoidectomy	ACT0106	1x6	20
	Haemorrhoidectomy	ACT0108	1x8	20
	Dialysis	ACT0303	3x3	5
	Vascular/Surgical	ACT0505	5x5	5
	Emergency/ Trauma	ACT0808	8x8	5
	Military	ACT1010	10×10	5



Anziam Bio Private Limited

 Bionest Bio-incubator, Regional Centre for Biotechnology, Gurgaon-Faridabad Expressway, Faridabad-121001
 Anziam Bio Pvt LtdE-127B, Ground FloorNear Vande Mataram Chowk, RIICO Industrial Area, Bhiwadi, District- Alwar, Rajasthan 301019

🔇 +91-7982045562

www.anziam.com

🔀 care@anziam.com

in linkedin.com/company/anziam-bio/

- facebook.com/Anziam-Bio/
- https://bit.ly/Anziam-Bio